


## Filling the medication cassette reservoir

This should be a clean/aseptic technique.

1. Open and inspect the cassette to ensure that it is intact and sterile.
2. Attach a Luer lock syringe containing the medication to the end of the line.
3. Tilt the cassette to approx 70° angle and fill. When 3/4 of the fluid is inserted, clamp the line. Gently tap the cassette to collect air in the outlet point.
4. Open the clamp and withdraw the air into the syringe.
5. Hold the syringe with the tip downwards and finish filling the reservoir.
6. Clamp the line and disconnect the syringe. Cap the line with a sterile cap, or attach an extension line.



## External reservoir adaptor

For infusion of larger volumes of fluid, an External reservoir adaptor (ERA) can be used to attach a non-vented fluid bag. This should be a clean/aseptic technique.




1. Open and inspect the ERA to ensure it is intact and sterile.
2. Close clamp to stop flow of fluid down the line.
3. Remove sterile cap and spike the fluid bag.

## Changing Medication Cassette / ERA

Follow your standard protocol for preparing the cassette.

1. Press . Pump displays **Stop pump?**
2. Press **Yes**. Pump displays **Pump is stopping**.
3. Clamp the tubing and disconnect from the patient.
4. Unlock the cassette latch and remove the empty cassette.
5. Attach a new cassette (see *Attaching cassette/ERA*)
6. **Set Reservoir Volume?** press **Yes**.
7. Enter security code (see Unlocking the pump)
8. Press to reset to the protocol volume displayed, or and to adjust.
9. When the reservoir volume is correct, press **Save**. Pump displays **Reservoir volume.....saving**.
10. Prime the tubing when prompted. (See Prime the tubing).
11. **Review pump setting and start pump?** Press **Yes**
12. Press **Review**, then check the patient specific parameters and Accept Value on each parameter.
13. When you have finished the review, press **Next** to continue. **Start pump?** Press **Yes**.



## Troubleshooting

### Air in-line. Press 'Acknowledge' then prime tubing

There may be air bubbles in the line; the pump will not run.

1. Press **Acknowledge** to silence the alarm.
2. Before clearing the air from the tubing, clamp the tubing and disconnect from the patient.
3. Remove the air using **Prime Tubing** on the **Tasks** menu, then reconnect to the patient.

### Battery low. Replace battery

The battery power is low but the pump is still operable.

1. Press **Acknowledge** to silence the alarm.
2. Change the batteries or recharge the battery pack as soon as possible.

### Downstream occlusion. Clear occlusion between pump and patient

High pressure may be due to blockage, kink in line or closed clamp.

- \* Remove the obstruction to resume operation without stopping the pump. or,
- \* Stop the pump, remove obstruction, review the program settings then restart the pump.

CADD®-Solis and the Smiths design mark are trademarks of the Smiths Medical family of companies. The symbol \* indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.

Smiths Medical MD, Inc.  
1265 Grey Fox Road, St. Paul, MN 55112, U.S.A.

European Representative:  
**Smiths Medical International Limited**  
Watford, Hertfordshire, WD24 4LG, U.K.  
www.smiths-medical.com  
Smiths Medical - part of the global technology business Smiths Group

**Smiths Medical International Limited**  
Military Road, Hythe, Kent, CT21 6JL, U.K.



© 2008 Smiths Medical family of companies.  
All rights reserved.

Part No. TPF-01132-101

smiths medical

smiths medical

## Quick Guide for the CADD®-Solis

Refer to the Operators Manual for the full instructions, including Warnings and Cautions, for the safe use of this product.



- Continuous infusion pain management with or without clinician bolus and/or patient controlled doses.
- Pump programmed by choosing a protocol: therapy, qualifier and drug, or programmed manually.

## Insert batteries

- Open the battery compartment and insert four AA batteries matching the + and - markings inside the battery compartment.



## Switch on

- To switch the pump on, press and hold the power switch.
- The pump carries out self-tests and sounds six beeps when the tests are complete.
- Screen shows **Start New Patient?**




- Close the compartment cover when the batteries are in place.

CADD

## Set up the pump for a new patient

Ensure you have the pump, a fresh set of 4 AA batteries, a filled medication cassette (or ERA and fluid bag). You also need to know the Keyboard Access code for the pump, or have the pump key.

1. Insert batteries and press the power switch to turn pump on.
2. **Do you want to start a new patient?** Press **Yes**.
3. Unlock the keypad.
4. Choose a **Therapy** and press **select**, **Qualifier** and press **select** and a **Drug** and press **select**.




5. **Are these correct?** Verify the chosen protocol and press **Yes**.
6. **Review pump settings.** Press **Review**.
7. Carefully check the patient specific parameters. Press **Accept Value** on each parameter, or press **select** to change.
8. When you have finished the review, press **Next** to continue.
9. **Cassette not attached. Attach cassette before starting pump.** Attach, latch and lock cassette to the pump.
10. **Prime Tubing?** Press **Prime**. Press **Stop Priming** when complete.
11. **Continue Priming?** Press **No**.
12. **Start Pump?** Press **Yes** when you have connected the tubing to the patient, and are ready to start the infusion.

## Change the program with the pump stopped



You can change settings when the pump is stopped. Before you restart the pump, the screen displays all the patient specific parameters which you must review and accept.

1. Use **▲** and **▼** to highlight the patient specific parameter you want to change, for example, Continuous Rate and press **select**.
2. Unlock the keypad.

The screen shows the current setting of your selected parameter. For rates, the range of soft limits for the protocol are shown in green, and the range of hard limits are in amber.




3. Use **▲** and **▼** to set the new value and press **Confirm**.
4. If a soft limit is exceeded, you see the message **Confirm soft limit override?** Press **Yes**.
5. When the change is complete, press **confirm**.
6. **Review pump settings.** Press **Review**.
7. Carefully check the patient specific parameters. Press **Accept Value** on each parameter, or press **select** to change.
8. When you have finished the review, press **Next**.
9. **Start Pump?** Press **Yes** to start the infusion. The Home screen is displayed, showing the **Reports** and **Tasks** menu keys.




## Attaching the cassette / ERA


Clamp tubing and check that the cassette lock is in the open position.



2. Insert the cassette hooks into the hinge pins on the bottom of the pump.




3. Push down on the latch as you move the cassette up into position. You feel a click as the cassette engages with the pump.
4. Lift the latch up into the closed position.
5. Insert the pump key and turn it clockwise to lock the cassette.



## Prime the tubing

Ensure that the pump is stopped, the tubing is disconnected from the patient, and the clamp is open on the patient tubing.


1. From the **Tasks** menu, use **▲** and **▼** to highlight the **Prime Tubing** option, then press **select**.
2. When prompted, unlock the pump.
3. Press **Prime** to start priming.
4. Press **Stop Priming** when the air is removed, or the delivery will stop automatically at 10 mL.



## Unlocking the pump

You can start and stop the pump, and deliver a PCA dose (if one is programmed) without unlocking the pump. For all other tasks, you must enter the security code. When the security code is required, the pump displays **Enter code to continue**.


1. Use **▲** and **▼** until the first digit of the code is shown on the screen.
2. Press **select** to advance to the next digit.
3. Repeat with the second and third digits, then press **select** or **Accept Value** to finish entering the code.





## Reports

The **Reports** function is on the Home screen and also appears as an option on the **Tasks** menu. Reports can be displayed when the pump is running, or when it is stopped.

1. Press **Reports** to display the **Select Report** menu.
2. Use **▲** and **▼** to highlight a report title and press **select**. The report is displayed on the screen.



Examples: PCA Dose Graph      Delivery History & Pie Chart



3. Press **Back** to return to the **Select Report** menu.