

Shiley Uncuffed Tracheostomy Tube

Tube Contains Metal

Great Ormond Street
Hospital for Children



NHS Foundation Trust

Made from opaque thermo sensitive PVC with a thin walled shaft, tapered tip and universal 15mm termination

SPECIAL INSTRUCTIONS

This tubes does **NOT** contain any metal and should be the tube used for MRI of the Head and Neck

Changed – PRN

Weekly Tube Changes

Shiley Tubes are not reusable and should **NOT** be sterilised and used again



Tube size fg **NEO/PAED** (delete as appropriate)

Suction Length cm Catheter Size fg

Last Tube Change .../.../...

Contact Jo Cooke ANP bleep 0712 or ENT On Call bleep 1020 if you have any concerns

STOMA CARE

Daily/PRN tape changes must be carried out.

Use of cotton tapes and Trachi-Dress

Correction tension – one finger space between tapes and child's neck

Emergency Box

Tube of the same size (Paed/Neo)

Tube one size smaller Suction Catheter (same size as suctioning)

KY Jelly

Tracheostomy Tapes

Round Ended Scissors